

### Silica Basin with Spout fused

Cat. No.	Cap ml
29	20
30	35
31	70
32	100
33	150
34	200
35	200

### Boats Combustion with Handle Fused silica


Cat. No.	Length	Width	Height
36	50	12	8
37	77	12	8
38	78	20	12
39	103	20	12


### Rubber Tubing For Laboratory Use

Cat. No.	Bore mm	W.T. mm
40	3.0	1.0
41	5.0	1.6
42	6.0	1.6
43	8.0	1.6
44	9.5	1.6
45	12.7	2.3
46	15.8	2.3
47	19.5	3.1
48	25.0	3.1

### Rubber Tubing For Pressure (VACUUM)

Cat. No.	Bore mm	W.T. mm
49	3.0	3.0
50	4.5	3.1
51	6.3	3.1
52	8.0	3.1
53	9.5	3.1
54	12.7	4.5
55	15.8	6.3
56	19.5	6.3
57	25.8	6.3





### Rubber Corks: acid and alkali proof with out holes

Cat. No.	Bottom Dia	Top Dia	Height
58	9	12.5	16 mm
59	12	16	18 mm
60	17	23	24 mm
61	21	27	28 mm
62	26	30	32 mm
63	32	43	41 mm
64	40	51	44 mm
65	51	60	42 mm
66	60	75	40 mm
67	73	80	40 mm

### Pippette Bulb Rubber

Cat. No.	Size
68	Small
69	Medium
70	Big

### Filteration Rubber Cone

Cat. No.	Size
71	Small
72	Medium
73	Big

### Filteration Rubber Bunks Solid

Cat. No.	Size. ml
74	250
75	500
76	1000
77	2000
78	5000
79	10000
80	20000

### Rubber Bellow Aquapet

Cat. No.	Size
81	Single Bellow
82	Double Bellow
83	Vaquapet small

### PTFE (Teflon) Tubing Lenght 3-8 Mtrs.

Cat. No.	ID mm	OD mm
84	1.6	3.2
85	4.8	6.4
86	6.4	8.0
87	8.0	9.5
88	11	12.7

### Transparent Silicone Rubber Tubing

Cat. No.	ID mm	OD mm
89	2	4
90	3	5
91	4	6
92	6	9
93	8	11
94	10	14
95	12	18
96	19.5	25.5

### Chemical Thermometers - Mercury in Classes Total Length 25 to 30 cms. Immersion 76mm - Zeal


Cat. No.	Ranges	Sub division
97	0 to 50 C	0.5 C
98	10 to 110C	1 C
99	0 to 150 C	1 C
100	0 to 250 C	1 C
101	0 to 360 C	2 C


### LABORATORY THERMOMETERS (Mercurial)

Cat. No.	Ranges	Length
104	0°-50°	30 cms
105	-10°-110°	30 cms
106	0°-150°	30 cms
107	0°-250°	30 cms
108	0°-360°	30 cms

### RED PETROLIUM THERMOMETERS (Alcohol)

Cat. No.	Ranges	Lengths
118	-50°+50°	30 cms
119	-50°+100°	30 cms
120	-100°+50°	30 cms
121	-100°+100°	30 cms





### McLeod Guage

Cat. No.	Ranges
123	10-0.01mm
124	20-1mm

Cat. No.	Ranges
125	50gms (6/Pack)

### Minimum immersion length (Long Stem)

Cat. No.	Ranges	Length
126	250	6" To 18"
127	110	6" To 18"
128	360	12" To 24"
129	250	12" To 24"
130	150	12" To 24"
131	110	12" To 24"

### Universal Clamps Cross Pattern Jaws Lined with Cork or Asbestos

Cat. No.	
132	M. S Casting Copper Oxidised
133	Alluminium Alloy Casting Powder Coated
134	Jumbo Size 60 -130 mm Powder Coated

### Three Finger Universal Clamps with Rubber Sleeve and Bossheads

Cat. No.	
135	M. S Casting Copper Oxidised
136	Heavy Bossheads 15 mm Hole 25 mm Length

### Boss Heads

Cat. No.	
137	M. S Casting Copper Oxidised 12 mm.
138	M. S Strip Type 12 mm
139	M. S Coasting Copper Oxidised 19 mm
140	M. S casting Copper Oxidised 25 mm
141	Boss Heads Special Types Block Model

### Retorte Stand with heavy M.S. Base and Rod

Cat. No.	Base x Rod Size mm
142	200 x 125 x 600
143	200 x 150 x 750
144	300 x 175 x 900
145	200 x 125 x 600 Powder coated
146	225 x 150 x 750 Powder coated
147	300 x 175 x 900 powder coated
148	200 x 125 x 600 S.S.304
149	225 x 150 x 750 S.S 304
150	300 x 175 x 900 S.S. 304



### Retorte Stand with heavy M.S. Base and Rod

Cat. No.	Base x Rod Size mm
151	V-Shape Heavy Duty 1.5 Mtr. Height
152	H-Shap Heavy Duty 1.5 Mtr. Height

### Retort Ring Fitted with Boss head (For Separating Funnels)

Cat. No.	Size mm
152/a	50
153	75
154	100
155	125
156	150

### Spatula S.S. Both End flat or one end Spoon

Cat. No.	Inches
157	4"
158	6"
159	8"
160	10"
161	12"

### Screw Clips Made of Thick Brass

Cat. No.	Size mm
162	12
163	18
164	25

### Laboratory Jackbox M.S. type


Cat. No.	Top Size mm
165	80 x 70
166	14 x 12
167	20 x 20
168	25 x 25
169	30 x 30

### Brushes with Nylon Bristles


Cat. No.	Type
170	Test Tube Brush
171	Bottle Brush
172	Pipette Brush
173	Burette Brush
174	Flask Brush

### Nickel Crucible with Lid Heaving 99.5% Purity Carbon Free Extra Thick

Cat. No.	MI.
175	15
176	25
177	50
178	100



For More Sizes See Next Page



**Nickel Crucible with Lid Heaving 99.5% Purity Carbon Free Extra Thick**

Cat. No.	MI.
177	50
178	100

**Nickel Spatula Rod Type**

Cat. No.	Inches
179	4"
180	6"
181	8"
182	10"
183	12"

**SS Foreceps**

Cat. No.	Inches
184	4"
185	6"
186	8"
187	10"
188	12"

**Nickel Basin with Spout Extra Thick**

Cat. No.	Dia mm
189	75

**Nickel Combustion Boat with Handle**

Cat. No.	Size mm
190	95 x 17 x 11

**Heating Mantles Single Size**

Cat. No.	Cap. ml
191	250
192	500
193	1 ltr
194	2 ltr
195	3 ltr
196	5 ltr
197	10 ltr
198	20 ltr
199	50 ltr

Note : All the above heating mantles are with energy Regulators

### Water Bath Double Walled Thermostatic Control Inside S.S. Vessel

Cat. No.	Size ml
200	250
201	500
202	1 ltr
203	2 ltr
204	3 ltr
205	5 ltr
206	10 ltr
207	20 ltr

### Oil Bath With Digital Temperature Controller

901	250
902	500
903	1 ltr
904	2 ltr
905	3 ltr
906	5 ltr
907	10 ltr
908	20 ltr

### Heating Coils SS Type

Cat. No.	Size ml
208	250
209	500
210	1 ltr
211	2 ltr
212	3 ltr
213	5 ltr
214	10 ltr
215	20 ltr

### Heating Coils Glass Type

Cat. No.	Size ml
216	250
217	500
218	1 ltr
219	2 ltr
220	3 ltr
221	5 ltr
222	10 ltr
223	20 ltr

### BEAKERS : Polypropylene

Cat. No.	Cap ml
224	50
225	100
226	250
227	500
228	1000
229	2000
229/A	5000



**WASH BOTTLES : Polyethylene**

Cat. No.	Cap ml
230	125
231	250
232	500
233	1000

**WASH BOTTLES(New Type) : Polyethylene**

Cat. No.	Cap ml
234	250
235	500
236	125

**FLASK STAND : Polyethylene**

Cat. No.	Particulars
237	Flask Stand

**DESICCATOR (vacuum) : Polyethylene**

Cat. No.	Size mm	Flange O.D. mm	I.D. mm
238	150	171	150
239	200	232	199
240	250	276	240

**Desiccator (Plain) : Polyethylene**

Cat. No.	Size mm	Flange O.D. mm	I.D. mm
241	150	171	150
242	200	232	200
243	250	275	240

**MEASURING CYLINDERS : Polypropylene**

Cat. No.	Cap ml
244	10
245	25
246	50
247	100
248	250
249	500
250	1000
251	2000

**MEASURING JUGS : Polypropylene**

Cat. No.	Cap ml
252	250
253	500
254	1000
255	2000



### FUNNELS - LONG STEM : Polyethylene

Cat. No.	Dia mm
256	50
257	65
258	75
259	100
260	150

### BUCHNER FUNNEL: Polyethylene

Cat. No.	Dia mm
261	70
262	110


### TEST TUBE STAND : Polypropylene

Cat. No.	Particulars
263	13 mm x 31 Tubes
264	16 mm x 31 Tubes
265	20 mm x 20 Tubes
266	25 mm x 12 Tubes
267	25 mm x 18 Tubes
268	32 mm x 12 Tubes
269	25 mm x 36 Tubes
270	13 mm x 18 Tubes
271	15 mm x 18 Tubes
272	18 mm x 12 Tubes
273	13 mm x 62 Tubes
274	16 mm x 62 Tubes
275	20 mm x 40 Tubes

### TEST TUBE STAND : Aluminum

Cat. No.	Particulars
276	13 mm x 31 Tubes
277	16 mm x 31 Tubes
278	20 mm x 20 Tubes
279	25 mm x 12 Tubes
280	13 mm x 62 Tubes
281	16 mm x 62 Tubes
282	20 mm x 40 Tubes
283	25 mm x 24 Tubes
284	25 mm x 18 Tubes
285	25 mm x 36 Tubes
286	32 mm x 12 Tubes





### Round Magnetic Stirring Bar with or without Pivot Ring

Cat. No.	Dia x Length
287	5 x 8
288	5 x 10
289	6 x 12
290	6 x 18
291	7 x 15
292	8 x 20
293	8 x 25
294	8 x 30
295	8 x 50
296	9 x 20
297	9 x 25
298	9 x 35
299	10 x 40
300	10 x 50
300/A	10 x 60
300/B	13 x 35
300/C	13 x 75

### Magnetic Stirring Bar Oval Shape Imported Model

Cat. No.	Dia x Length
301	4 x 8
302	5 x 10
303	6 x 15
304	10 x 20
305	12 x 25
306	16 x 30
307	16 x 35
308	20 x 40
309	20 x 50
310	20 x 64
311	20 x 70

### Octagonal Magnetic Stirring Bar

Cat. No.	Dia x Length
312	8 x 13
313	8 x 16
314	8 x 25
315	8 x 38
316	9 x 13
317	9 x 19
318	9 x 25
319	9 x 38
320	9 x 51
321	13 x 15
322	13 x 75

### Magnetic Stirring Bar Egg Shape Indian Model

Cat. No.	Dia x Length
323	5 x 10
324	6 x 15
325	10 x 20
326	12 x 25
327	16 x 30
328	16 x 35
329	20 x 40
330	20 x 50
331	20 x 70

### Lab Bottles Screw Type MURTHY'S Complete Set.

Cat. No.	Size ml.
339	100
340	250
341	500
342	1000
343	2000
344	5000
345	10000


Cat. No.	
346	Alluminium Cap
347	PP Cap
348	4 Holes Rubber Cork
349	2 Holes Rubber Cork
350	Plain Rubber Septum
351	Teflon Coated Septum
352	PP Caps Turning only

### Turn on Rubber Septums Latex White


Cat. No.	Cone Size
353	14/23
354	19/26
355	24/29
356	29/32
357	34/36

### Turn on Rubber Septums silicone Red

Cat. No.	Cone Size
358	14/23
359	19/26
360	24/29
361	29/32
362	34/36



Note : Set Contains of Alluminium Cap, PP Cap, 4 Holes Rubber Corck, Plain Rubber Septum. On enquiry : 2 Holes Rubber Corck, Teflon Coated Rubber Septum.



**Rotavapor condenser (For Buchi & Heidolph Models)**

Cat. No.	
363	condenser
364	condenser vacuum key
365	condenser socket & key

**Homogenizer with Teflon Handle**

Cat. No.	Cap ml
366	2
367	5
368	15
369	25
370	50

**Liophylizer Tube one mark 50mm from bottom**

Cat. No.	Dia x Length
371	9 x 170

**RB Base Made of Silicone Rubber**

Cat. No.	Cap ml
393	100
394	250
395	500
396	1000
397	2000
398	3000
399	5000
400	10 ltr
401	20 ltr

**Vials screw type (Ambered) Neutral Glass**

Cat. No.	Cap ml
402	2
403	5
404	10
405	15

**Vials screw type (Clear)**

Cat. No.	Cap ml
406	2
407	5
408	10
409	15

**Vials screw type (Amber Clear Soda Glass)**

Cat. No.	Cap ml
410	60
411	100

For Borosile glass please send your enquires.